# IKA

## Screening system

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**About this instruction manual**

- ▶ Before commissioning, read and observe the safety instructions.
- ▶ Keep the operating instruction manual accessible for all users.
- ▶ Read and observe the operating instructions for the power supply unit the Magnetic Stirrer RCT digital.

**Accompanying documents**

Operating instructions for:

- Magnetic Stirrer RCT Digital
- Rohde&Schwarz HMP4040 power supply unit

**Intended purpose**

The IKA Screening-systems, 6 cells or 8 cells, are used to perform electrosynthesis for research and development or small production runs in laboratories. The multiple-cell-design allows for testing multiple conditions in the same vessel holding different parameters constant.

Safe functioning is only assured with the accessories listed in this instruction manual or recommended by the manufacturer.

Which aids and operating material are used is the responsibility of the user.

Damaged or incorrectly fitted parts can cause serious accidents. Only use devices in full working order. Only use devices in their original state.

Devices must not be opened unless instructed to do so in the manual of said device.

The cells in the Screening system must not be operated in explosive atmospheres or with materials that could form an explosive mixture.

**Other risks****Danger of injury due to electric shock from the power supply unit**

- ▶ Do not use the power supply unit if it is not working correctly. The device guard may be damaged. If you are unsure whether the power supply unit is working correctly, have it repaired by a specialist workshop.
- ▶ Check supply lines for damage to insulation or exposed metal before commissioning. Replace damaged supply lines with new ones.
- ▶ Exceeding the Low Voltage Protection!  
For the series connection of all output voltages, it is possible to exceed the low voltage protection of 42 V. Please note that in this case any contact with live components is life-threatening. It is assumed that only qualified and trained personnel service the power supplies and the connected loads. When changing the mains voltage, it is essential to replace the fuse. Otherwise the device may be destroyed.

**Danger of explosion due to sparks from the power supply unit**

- ▶ Do not use the device in hazardous areas, it is not EX-proof.
- ▶ In the case of substances, which can form an ignitable mixture, appropriate protective measures, such as working under a fume hood, must be taken..

**Danger of injury due to tipping or falling parts**

- ▶ Ensure that the parts are set up so they cannot tip over start to move.
- ▶ Do not exceed permitted loads, torques and separations.
- ▶ Before every use check whether all (clamp) connections are firm and secure.
- ▶ Check that the caps of the cells and electrodes are properly tightened before each use.
- ▶ The laboratory device may only be moved on the stand when the device is at a standstill and unplugged.

**Danger of injury due to inhalation of substances damaging to health**

Contact with lead, nickel, cobalt, and lead-bronze electrodes can lead to poisoning.

In order to avoid injury to persons and damage to property, please observe the relevant protective and accident prevention measures when processing hazardous substances.

**Requirements for operators / users**

- ▶ Ensure that only authorised and trained persons work with the devices.
- ▶ Wear personal safety equipment in accordance with the hazard class of the devices.
- ▶ Observe GLP guidelines (Good Laboratory Practice).
- ▶ Repairs must only be carried out by a qualified technician who is aware of the possible dangers.
- ▶ Place devices on an even, stable, clean, non-slip, dry, fire-resistant surface.
- ▶ Before setting up devices ensure that there is 1 m<sup>2</sup> of free space.
- ▶ Pack devices carefully and check for damage. In case of damage in transit, you should notify the carrier immediately.

**Requirements for use in Japan**

- ▶ When connecting the device to the power, a voltage convertor with step-up is required. Both power supplies can be plugged into the provided voltage convertor.
- ▶ Make sure to switch to the correct voltage at the device power supply unit.
- ▶ Please refer to the operating manual of the power supply unit.

## Symbols used in this instruction manual

“Safety” chapter

“Scope of delivery” chapter

“Assembly” chapter

“Operation” chapter

“Fault rectification” chapter

“Cleaning and maintenance” chapter

“Technical data” chapter

### Item numbers

Indicates relevant device components for actions.

### Correct / result

Indicates correct implementation or the result of an action step.

### Incorrect

Shows incorrect implementation of an action step.

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See more accessories on [www.ika.com](http://www.ika.com).

<p>Geteilte Zellen / Divided Cells</p>	
<p>Ungeteilte Zellen / Single Cells</p>	
<p>Aluminium Block (6 Zellen) / Aluminum block (6 cells)</p>	
<p>Aluminium Block (8 Zellen) / Aluminum block (8 cells)</p>	
<p>Kronenrührer / Crown stirrer</p>	
<p>Elektroden / electrodes</p>	
<p>RCT Digital</p>	
<p>Netzgerät / Power supply</p>	
<p>Stativstab H 16V / Support rod H 16V</p>	
<p>Haltestange H38 / Holding rod H38</p>	
<p>Kreuzmuffe H44 / Cross sleeve</p>	

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
1. Setting up and connecting the magnetic stirrer
2. Setting up and connecting the power supplies
3. Assembly of the reaction chambers:
  - a. split chamber: check tightness of the screws holding the chambers together
  - b. Add stirrer with the crown pointing downwards
  - c. Mount electrodes in cover, dipped sufficiently, tighten screws by hand
  - d. Mount cover with electrode
4. Place the reaction blocks on the magnetic stirrer
5. Fill the reaction chambers outside or inside the reaction block, pay attention to the level. Carefully place the electrodes
6. Connect the cables from power supply to electrode  
Bottom row  
(red/red)  
Channel 1 -> reaction chamber 1  
Channel 2 -> reaction chamber 2  
Etc.
7. Set RCT temperature and speed
8. Remote control via labworldsoft (serial cable to RCT / power supply) or via USB

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**Note!**

Read the mounting and safety instructions of the **IKA** boss head clamp prior to using it.

**The device must not be suspended from the support rod!**

**Note!**

When using attachments with a diameter over 180 mm, use support rod in conjunction with an extension. (see "Accessories")

**Risk of tipping!**

Make sure that the center of mass of the attached device does not protrude over the safe area indicated below with a dotted rectangle.

**Procedure**

1. Remove the threaded plug in the threaded hole.
2. Remove the protective cap from the support rod.
3. Place the washer between the housing and the nut.
4. Screw in the support rod manually until it no longer can be tightened.
5. Tighten the nut using an SW17 open-end spanner.
6. Assemble the accessories using a boss head clamp.

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- ▶ Tighten screws.
- ▶ Check seals.

### Power supply

For a complete fault correction list, see the manufacturer's operating instructions.

- ▶ Ensure that the fault is related to the device itself and that external connections are not faulty.
- ▶ Ensure that the device has been precisely calibrated during the last year.

If the device is not in working order, carry out the following checks:

- ▶ Ensure that the AC power cable is connected to the power supply unit.
- ▶ Ensure that the mains switch on the front panel is pressed down.
- ▶ Ensure that the correct mains fuse has been fitted.
- ▶ Ensure that the electrodes are correctly connected to the clamps.

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### Electrodes

- ▶ Visual inspection
  - ▶ In case of deposits clean by hand
- Appropriate cleaning agents: observe material combinations (harmless, e.g. water, acetone, acetonitrile)

### Power supply

For a complete list of maintenance and cleaning regulations, see the manufacturer's original operating instructions.

To avoid electric shocks, disconnect the device from the AC power mains power supply and remove all measuring lines before cleaning. Clean the surfaces of the device with a soft, lint-free, slightly moist cloth.

- Do not use any cleaning agents or solvents.
- Never clean inside.

### Disposal

Devices must be fully disassembled and the materials removed and disposed of separately in accordance with applicable regulations (e.g. separate disposal of oils, plastics, metal and electronic industrial waste).

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## Technical data

### Heating block for single cells IKA-Item: 20015943

#### Chamber for cells

Number of chambers for cells		8
Chamber diameter	mm	25
Chamber depth	mm	30
Diameter	mm	150
Material		Aluminum, anodized
Weight	kg	1.4
Max. temperature	°C	200

#### Hole for temperature probe

Diameter	mm	3.3
Depth	mm	30

#### Cell

Diameter	mm	25
Material		PTFE
Weight	kg	0.1
Volume	ml	20

#### Stir bar

Type		Single cross head
Diameter	mm	14
Height	mm	12
Material		PTFE

### Heating block for divided cells IKA-Item: 40001854

Diameter	mm	150
Material		Aluminum, anodized
Weight	kg	2.1
Max. temperature	°C	200

#### Chamber for cells

Number of chambers for cells		6
Chamber diameter	mm	50
Chamber depth	mm	60

#### Hole for temperature probe

Diameter	mm	3.3
Depth	mm	30

#### Cell

Diameter	mm	50
Material		PTFE
Weight	kg	0.2
Volume	ml	2 x 10

#### Frit

Material		Glass
Diameter	mm	10
Porosity	mm	14

#### Stir bar

Type		Single cross head
Diameter	mm	14
Height	mm	12
Material		PTFE

**IKA-Werke GmbH & Co.KG**

Janke & Kunkel-Str. 10

D-79219 Staufen

Tel. +49 7633 831-0

Fax +49 7633 831-98

sales@ika.de

**www.ika.com**

40004123\_20000021616\_OMI\_screening\_system\_DE-EN\_1.1.docx