

# Simple and Affordable



### One Titrator for Many Needs

This versatile titrator is ready to carry out all your acid/base inflection point, pH endpoint or precipitation/redox indicative titrations. Be prepared for all basic titrations whatever your industry.



### Easy App Guidance

This titrator has been especially designed for simplicity. New operators are already familiar with the Smartphone like Apps oriented user interface. Start your routine analysis with a LongClick™ on the home screen. The intuitive navigation speaks your language and 14 others. Learning and using a new analytical instrument has never been so easy.



### Swiss Engineered

EasyPlus™ Titrators are designed and constructed by our Swiss development team according to the highest possible software and hardware standards. Their renowned expertise has made the impossible possible; to develop such a high quality instrument with such a small footprint.



### iTitrate for Easy Programming

iTitrate Intelligence makes method programming simple with EasySetup guidance. Set only the minimal parameters and the titrator programs the rest, including correct calculations based on desired measurement units. Easy setup equals swift implementation.



### Easy Pro Titrator

The all in one Titrator

This versatile titrator is suitable for a range of determinations making it ideal for simple titrations of most sample types. The simple 'App' style user interface reduces training time while increasing work speed and simplicity. iTitrate Intelligence's guided set-up assists in method programming to ensure samples are tested consistently, regardless of operator, and with minimal effort. Eliminate manual calculations and documentation with a printer or EasyDirect™ Titration PC Software.

The EasyPlus Titrator is capable of performing titrations determined by:

- pH endpoint titration in juice, wine, p- and m-value
- Acid/base inflection/equivalence point titrations for free fatty acids in food and alkalinity etc..
- Precipitation titrations, such as chloride in food or water and sulfides in oils.
- Redox titrations, such as Vitamin C in beverages, and reducing sugar in wines, sulfite in water, etc..

# Easy Pro

## Standard delivery

Product	Quantity	Order Number
EG11-BNC Sensor pH aqueous	1	30043103
EasyPlus burette 20 mL	1	30043901
EasyPlus tubing set instrument	1	30065464
EasyStir GT	1	30065467
EasyPlus titration head GT	1	30041102
EasyPlus drying tube	1	30044701
Valve	1	30042860
Bottle head	1	30060023
Insert set titration head GT	1	30057635
Burette cover set	1	30057633

## Accessories

Product	Order Number
EasyPlus burette 10mL	30043900
EasyPlus burette 50mL	30043902
EG13-BNC sensor pH non-aqueous	30043104
EM40-BNC sensor redox mV	30043106
EM43-BNC sensor redox Ipot	30043105
EM45-BNC sensor precipitation	30043107
Bottle adaptor China	30060026
Bottle adaptor India	30060024
Bottle adapter Merck	23774
Bottle adaptor Fischer	23787
EasyPump	30064218
Siphon tips (5 pcs.)	23240
Glass beakers (80 mL), 20 pcs.	101446
Plastic beaker (80 mL), 120 pcs.	51109388
Temperature sensor 30K	51300164
EasyDirect instrument license for 3 titrators per PC	30065449
USB-P25 Strip Printer	11124301

## Technical specifications

pH/mV sensor input	Measuring range	± 2000 mV / 0 – 14 pH
	Resolution	± 0.1 mV
Polarized sensor Input	Current source redox Ipot	0 – 5 µA
	Resolution	± 0.1 µA
	Measuring range (voltage input)	0 – 2000 mV
Temperature sensor input	Measuring range	0 – 100 °C
	Resolution	± 0.1 °C
Power supply	Input voltage	24 VDC / 1.25 A
	Line input voltage	100 – 240 V ±10 %
	Input frequency	50/60 Hz
Environment	Ambient temperature	5 °C – 40 °C
	Atmospheric humidity	Max. 80% at 31 °C
Titrator dimensions	Width x depth x height / weight	170 x 220 x 350 mm / 1.9 kg (Plus KF stand 0.8 kg)
Control and display	Touchscreen	4.3 inch, 480 x 272 pixel RGB
Communication	USB-A	USB full / low speed
	USB-B	USB full / low speed
Application*	Repeatability	0.5% RDS
	Titration type	Direct, back compensated, blank compensated

\*Depends on the environment and handling

## Support or Service? No problem



### Application know-how

Would you like application know-how to titrate a specific sample?



### Expertise & help

Would you like best practice titration tips or require technical support expertise?



### Maintenance & verification

Do you want learn more about preventive maintenance or need to carry out a performance verification of your titrator?

Simply register to the EasyPlus™ Titration Webpage and quickly find comprehensive help for your daily tasks.

► [www.mt.com/easyplustitration](http://www.mt.com/easyplustitration)



Quality certificate ISO9001  
Environment certificate ISO14001  
Internet: <http://www.mt.com>  
Worldwide service

Subject to technical changes  
© 01/2013 Mettler-Toledo AG  
MarCom Analytical, 30073726

[www.mt.com/easyplustitration](http://www.mt.com/easyplustitration)

For more information