

# SCINTECK INSTRUMENTS LLC, USA

We carry complete inventory of Deuterium/Hollow Cathode and Tungsten Halogen lamp for all Equipment Manufactures  
SCINTECK Manufactures and supply replacement deuterium lamps for virtually all instrument manufacturers,


ABI, AGILENT, ANALYTIK JENA, BECKMAN, BEIJING PURKINJE, BEIJING BEIFEN-RAILI, BIO-RAD ERMA, HACH, HITACHI, JASCO, KNAUER, LINEUR, PERKIN ELMER, PG INSTRUMENTS, SCINCO, SCIEX, SHIMADZU, SHANGHAI WUFENG, SHANGHAI METASH, SPECTRA PHYSICS, SSI, TECHCOMP, THERMO, DIONEX, VARIAN, WATERS AND MORE

Sooner or later, every deuterium /Hollow Cathode/Tungsten Halogen reaches the end of its useful life. Many instrument manufacturers demand extremely high prices from their users for replacement lamps. However, as only very few of these companies produce deuterium lamps themselves, they generally outsource them from OEM Lamp manufacturers:

As a Manufacturer and supplier of these products and a close partner of other companies, who fit Hamamatsu lamps with instrument-specific brackets and flanges, we can offer prealigned lamps for almost all instruments at prices well below those demanded by the instrument manufacturers. We offer our customers competent technical advice and maintain stock of many common lamps, enabling us to deliver fast.

Our databases are compiled with great care and attention. We guaranty that the listed lamps are suitable for your instruments. Should this not be the case, we will take them back and refund the full purchase price. You, our customers, are working with cutting-edge technology, so please: if you are using instrumentation that isn't listed on our website or if you need a lamp that isn't shown there, please contact us with as much information as possible. We will do our best to help you.

Deuterium lamps emit an almost continuous spectrum of light ranging from the main UV wavelengths of 160 - 400 nm to the visible spectral range (800 nm). This makes them the ideal light source for high precision absorption measurements, e.g., in HPLC. SCINTECK offers a wide range of deuterium lamps: Lamps with a guaranteed lifetime of 2000 hours, with quartz or UV-glass envelopes, shine-through lamps, and different sorts of special lamps.



Scintek Instruments LLC, USA  
(571) 426-3598, (571) 426-1023, [sales@scintek.com](mailto:sales@scintek.com)

*scintek.com*